

emergency



# MEDUMAT Easy CPR

*The Guiding Emergency Ventilator*

# MEDUMAT Easy CPR

*The Guiding Emergency Ventilator*



## Safety is on Duty

Although patient resuscitation is not part of your daily work, you are regularly confronted with this challenging situation. Then it is all the more important to have a partner on your team that's completely reliable and capable of helping you with ventilation and Cardiopulmonary Resuscitation (CPR).

MEDUMAT Easy CPR was specially developed for emergency ventilation during and after resuscitation. You can, for example, trigger a ventilator breath right from the mask.

The ventilator ensures the highest degree of safety by giving you clear voice prompts that lead you through CPR procedures.

### **Intuitive operation, small and lightweight**

In developing MEDUMAT Easy CPR, we concentrated on ergonomic design, intuitive operation and a size and weight that would be easy to transport.



### fast resuscitation

- quick setting
- metronome > 100/min
- ventilation in 5 seconds

### ergonomic

- simple operation
- low weight
- easy to transport

### safe resuscitation

- no hyperventilation
- low risk of stomach inflation

## Increased safety and better outcomes

You can quickly make the main settings (frequency, tidal volume, pressure limits) which comply with the ERC 2010 Guidelines.

With the help of the integrated metronome and the voice prompts, you will find it easy to provide the recommended rate of 110 chest compressions per minute and two breaths within five seconds.

With the administration of a constant volume and a regular ventilation frequency, you can also avoid hyperventilation of the patient. The adjustable pressure limits help you to prevent insufflation of the patient's stomach.

**"The [...] Medumat Easy volume-cycled ventilators were technically efficient and reliable; their gas consumption was less than that of conventional emergency ventilators."\***

\*L'Her, E., Roy, A.: Bench Tests of Simple, Handy Ventilators for Pandemics: Performance, Autonomy and Ergonomy, Respiratory Care, Volume 56, Number 6, June 2011



# MEDUMAT Easy CPR – Cleverly Designed

## ERC 2010 guideline compliance

- Tidal volume:
  - 6–7 ml/kg weight
  - 500 ml – 600 ml at  $f = 10/\text{min}$
- Acoustic metronome signal  $f > 100/\text{min}$  ensures the correct chest compression frequency
- Ventilation phase in CPR mode: max. 5 seconds
- Both hands free for mask ventilation

## User-friendly operation

- Small and space-saving
- Light and handy
- Easy to transport
- Simple integration in CPR process

EN 1789

## Flexibility

- You decide in CPR mode when the patient is to be ventilated.
- Intuitive operation minimizes training needs.
- Portable system LIFE-BASE makes MEDUMAT Easy CPR mobile.
- Expandable functionality with Weinmann modules
- Practical accessories, can be safely stored

## Safety

- Clear and explicit voice prompts
- Metronome function ensures correct frequency for CPR
- Maximum ventilation pressure can be set
- Switch for ventilation pressure limit (20 mbar/45 mbar)
- Visual and acoustic alarms
- Safe selection of ventilation parameters
  - Dial for tidal volume and ventilation frequency



# LIFE-BASE

## MEDUMAT Easy CPR Designed for Mobility and Transport

MEDUMAT Easy CPR is available on many LIFE-BASE portable systems and offers you the following benefits in ventilation and oxygen therapy:

Easy to transport

Expandable functionality with function modules

Safe and practical storage of accessories



MEDUMAT Easy CPR on LIFE-BASE light (WM 28235)



MEDUMAT Easy CPR on LIFE-BASE III (WM 9155)



LIFE-BASE III MANV with MEDUMAT Easy CPR, MODUL Oxygen and 3-way distributor (WM 8000)

# MEDUMAT Easy CPR Speaks your Language

Professionals around the world depend on the helpful voice prompts from MEDUMAT Easy CPR, which is available with more than 20 languages:

### MEDUMAT Easy CPR, WM 28140

- |            |           |
|------------|-----------|
| German     | Danish    |
| English    | Norwegian |
| French     | Finnish   |
| Italian    | Icelandic |
| Spanish    |           |
| Portuguese |           |
| Dutch      |           |
| Swedish    |           |

### MEDUMAT Easy CPR, WM 28160

- |          |            |
|----------|------------|
| German   | Turkish    |
| English  | Indonesian |
| Russian  | Thai       |
| Polish   | Farsi      |
| Czech    | Persian    |
| Chinese  | Korean     |
| Japanese | Hindi      |
| Arabic   |            |



# This is How **MEDUMAT Easy CPR** Helps You – A Case Study

Just before the end of the working day, Hans S. doesn't feel well. But there's another pallet to move to the warehouse. Suddenly his heart stops beating. His breathing ceases. Everything goes dark.



The First Responders reach Hans S.

Now MEDUMAT Easy CPR helps to get ventilation and resuscitation started at once. It tells the rescuers what to do:



Voice prompt from MEDUMAT Easy CPR:

**“Adjust settings and connect patient”**



**“Perform chest compressions now”**

The metronome helps by giving a frequency of > 100/min.



**“Provide two ventilations now”**

Ventilation can be manually triggered. The First Responder has both hands free to position the mask.

## Technical data MEDUMAT Easy CPR



<b>Product class as per directive 93/42/EEC:</b>	IIb
<b>Protection from splashing water:</b>	IP54
<b>Standards used:</b>	EN 794-3, EN 60601-1, EN 1789
<b>Dimensions (W x H x D):</b>	100 x 145 x 90 mm incl. connections
<b>Weight:</b>	about 0.6 kg
<b>Temperature</b>	
■ for operation:	-18 °C to +60 °C
■ for storage:	-40 °C to +70 °C
<b>Gas supply:</b>	medical grade oxygen
<b>Gas consumption for internal management</b>	0 l/min
<b>Operating pressure:</b>	2.7 to 6 bar*
<b>Power supply</b>	
■ Lithium battery (accessible from outside device):	3.6 V; 5.2 Ah, expected service life > 2 years
<b>Setting options:</b>	- Deactivation of voice prompts - Deactivation of CPR metronome - Depending on device version, more than 12 different languages

<b>IPPV mode</b>																			
■ Ventilation frequency:	infinitely adjustable from 10 to 25 min <sup>-1</sup>																		
■ Tidal volume (V <sub>t</sub> ):	infinitely adjustable from 65 to 950 ml l																		
■ Correlation of tidal volume to ventilation frequency:	<table border="1"> <tr> <td>f (min<sup>-1</sup>)</td> <td>25</td> <td>20</td> <td>15</td> <td>12</td> <td>10</td> <td>10</td> <td>10</td> <td>10</td> </tr> <tr> <td>V<sub>t</sub> (ml)</td> <td>65</td> <td>100</td> <td>150</td> <td>300</td> <td>500</td> <td>600</td> <td>800</td> <td>950</td> </tr> </table>	f (min <sup>-1</sup> )	25	20	15	12	10	10	10	10	V <sub>t</sub> (ml)	65	100	150	300	500	600	800	950
f (min <sup>-1</sup> )	25	20	15	12	10	10	10	10											
V <sub>t</sub> (ml)	65	100	150	300	500	600	800	950											
■ Max. ventilation pressure:	20 or 45 mbar*																		
■ Double stroke process:	upon reaching max. ventilation pressure (e. g. with chest compressions)																		
<b>Demand flow mode</b>																			
■ Trigger:	< 1 mbar*																		
■ Peak flow:	47 to 52 l/min																		
■ Shut down pressure:	3 mbar*																		
<b>Visual alarm with voice prompt and alarm signal in event of:</b>	- Stenosis - Disconnection/ Apnoe - Decrease in supply pressure (< 2,7 bar) - Battery discharged/Device failure																		

\* 1 mbar  $\Delta$  1 hPa

## Accessories and replacement parts

1 EasyLung-Test lung	WM 28625	■ Lip membrane Set I (25 x WM 3211)	WM 15320
2 PEEP valve with connection nozzle, infinitely adjustable from 2.5 to 10 mbar (22 mm external diameter)	WM 3215	■ Lip membrane Set II (50 x WM 3211)	WM 15330
3 Patient hose system with patient valve for MEDUMAT Easy CPR	WM 22520	■ Membrane for expiration tube on patient valve	WM 3212
■ Ventilation tube, double lumen for MEDUMAT Easy CPR	WM 22647	■ Insert for spontaneous breathing tube	WM 3282
4 Patient hose system with patient valve, disposable, for MEDUMAT Easy CPR	WM 28110	■ Membrane for spontaneous breathing tube, red	WM 3284
■ Disposable ventilation hose set 10 x WM 28110	WM 15454	5 Reusable protective sleeve for ventilation hose for MEDUMAT Easy CPR	WM 8297
■ Disposable ventilation hose set 25 x WM 28110	WM 15455	6 Retrofitting set for device outlet filter for MEDUMAT Easy CPR	WM 15780
■ Disposable ventilation hose set 50 x WM 28110	WM 15456	■ MEDUtrigger for manual breath initiation for MEDUMAT Easy CPR	WM 20900
■ Patient valve for MEDUMAT Easy CPR, consisting of: Patient hose connection (WM 3213), Upper control element (WM 3181), Lower control element with spontaneous breathing option (WM 3285), Lip membrane (WM 3211)	WM 3280	■ Protective cap MEDUMAT Easy CPR device socket for MEDUtrigger	WM 28144
■ Lip membrane	WM 3211	■ Bacteria filter	WM 22162
		■ Test set for MEDUMAT Easy CPR	WM 15357
		■ Battery for MEDUMAT Easy CPR	WM 28045
		■ Filter cap for dust filter	WM 4954
		■ Dust filter for all MEDUMAT devices	WM 7571



## Weinmann medical technology

### Simply professional

You can rely on Weinmann. We set the standards for mobile solutions in patient ventilation. Our 130 years of experience and close association with professional users in emergency medical services, armies/public authorities and hospitals have led to the development of intelligent and perfectly coordinated technologies. In conjunction with the professionals in emergency, transport and disaster medicine, we pursue our shared goal of saving human lives.

**Would you like to know more?  
Call us on  
+49-(0)40-5 47 02-170**

### Germany

Weinmann Geräte für Medizin GmbH+Co.KG  
PO Box 540268 ■ 22502 Hamburg  
Kronsaalsweg 40 ■ 22525 Hamburg  
weinmann.de

T: +49-(0)40-5 47 02-0 Reception  
F: +49-(0)40-5 47 02-461 Reception  
T: +49-(0)40-5 47 02-120 Customer Service

Center for Production, Logistics, Service  
Weinmann Geräte für Medizin GmbH+Co.KG  
Siebenstücken 14 ■ 24558 Henstedt-Ulzburg

### France

Weinmann S.A.S. – Paris-Igny  
T: +33-(0)1 69 35 53 20  
weinmann-france.fr

### Russia

Weinmann SPb GmbH – St. Petersburg  
T: +7-(812)6-33 30 82  
weinmann.ru

### China

Weinmann GmbH+Co.KG (Shanghai Rep. Office)  
T: +86 21 32 26 26 98  
weinmann.cn

### Singapore

Weinmann GmbH+Co.KG (Singapore Rep. Office)  
T: +65 650 944 30  
weinmann.de

### Venezuela

Weinmann GmbH+Co.KG (Latin America Rep. Office)  
T: +58 274 416 44 66  
weinmann.de

